



Confirmation

Herewith i confirm, that the Erbe smoke evacuation system IES2 (10322-000) is the only validated smoke evacuation system in combination with the Erbe Electrosurgery systems VIO 3 (10160-000), VIO 300D (10140-100) and VIO 300S (10140-300) over the CAN-Bus connection (Erbe communication bus)

This enables:

- the activation of mode specific settings of the IES2
- socket specific settings of the IES2
- individual basic suction of the IES2
- individual run-on time of the IES2



Erbe Elektromedizin GmbH
Walchoernlestrasse 17
72072 Tuebingen, Germany

Michael Reick
Director Business Management